

## Cardiac Troponin I antibody [C5]

Cat. No. GTX10233

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, IP, Immunoassay, IHC, Purification
Reactivity	Human, Mouse, Rat, Rabbit, Goat, Bovine, Cat, Dog, Pig, Fish

Package  
200 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IP	Assay dependent
Immunoassay	Assay dependent
IHC	Assay dependent
Purification	Assay dependent

Not tested in other applications.

Calculated MW 24 kDa. ( [Note](#) )

Product Note Cross-reactivity with skeletal Tnl is more than 50%.

## Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Free human cardiac troponin I, native human cardiac troponin complex, or human cardiac troponin I peptides.
Purification	Protein A purified
Conjugation	Unconjugated

For full product information, images and publications, please visit our [website](#).

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).