

## Cardiac Troponin I antibody [C5]

**Cat. No. GTX10233**

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG2b
<b>Applications</b>	WB, IP, Immunoassay, IHC, Purification
<b>Reactivity</b>	Human, Mouse, Rat, Rabbit, Goat, Bovine, Cat, Dog, Pig, Fish

Package  
200 µg

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IP	Assay dependent
Immunoassay	Assay dependent
IHC	Assay dependent
Purification	Assay dependent

Not tested in other applications.

<b>Calculated MW</b>	24 kDa. ( <a href="#">Note</a> )
<b>Product Note</b>	Cross-reactivity with skeletal Tnl is more than 50%.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.09% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
<b>Concentration</b>	Batch dependent (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Free human cardiac troponin I, native human cardiac troponin complex, or human cardiac troponin I peptides.
<b>Purification</b>	Protein A purified
<b>Conjugation</b>	Unconjugated



For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

**Note**

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).