

Cardiac Troponin I antibody [M155]

Cat. No. GTX10237

| | |
|-------------|--|
| Host | Mouse |
| Clonality | Monoclonal |
| Isotype | IgG1 |
| Application | WB, IP, Immunoassay, IHC, Purification |
| Reactivity | Human |

Package
200 µg

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB | Assay dependent |
| IP | Assay dependent |
| Immunoassay | Assay dependent |
| IHC | Assay dependent |
| Purification | Assay dependent |

Not tested in other applications.

Calculated MW 24 kDa. ([Note](#))

Product Note No cross-reactivity with skeletal muscle Tnl.

PROPERTIES

| | |
|---------------|---|
| Form | Liquid |
| Buffer | PBS |
| Preservative | 0.09% Sodium azide |
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.) |
| Immunogen | Free human cardiac troponin I, native human cardiac troponin complex, or human cardiac troponin I peptides. |
| Purification | Protein A purified |
| Conjugation | Unconjugated |

For full product information, images and publications, please visit our [website](#).

Date 2024 / 05 / 17 Page 1 of 2

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).