

Albumin antibody [1A9]

Cat. No. GTX10242

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ELISA, Sandwich ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent
Sandwich ELISA	Assay dependent
Note : Capture : GTX10241/GTX10242, Detection : GTX10242/GTX10241	

Not tested in other applications.

Calculated MW	69 kDa. (Note)
Product Note	Does not react with bovine serum albumin, egg white albumin or human alpha fetoprotein.

Properties


Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Human serum albumin
Purification	Protein A purified
Conjugation	Unconjugated

For full product information, images and publications, please visit our [website](#).


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTXT10242 ELISA Image

Titration of human serum in two albumin sandwich immunoassays. Human Albumin antibody [15C7] (GTXT10241) for capture and Human Albumin antibody [1A9] (GTXT10242) for detection


GTXT10242 ELISA Image

Titration of human serum in two albumin sandwich immunoassays. Human Albumin antibody [15C7] (GTXT10241) for capture and Human Albumin antibody [1A9] (GTXT10242) for detection



For full product information, images and publications, please visit our [website](https://www.genetex.com).