

ERAB antibody [10D12]

Cat. No. GTX10259

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, IHC-Fr, Dot
Reactivity	Human

References (1)
Package
100 µl

Applications

Application Note

Dot: Use at a dilution of 1/100 - 1/1000. IHC: Use at a dilution of 1/100 - 1/1000. WB: Use at a dilution of 1/100 - 1/1000. Predicted molecular weight: 27 kDa. Optimal dilutions/concentrations should be determined by the end user.

Calculated MW 27 kDa. ([Note](#))

Product Note This antibody recognizes the recombinant human ERAB protein (27kDa) and the endogenous ERAB protein in cell extracts in immunoblots and ELISA assays, with virtually no crossreactivity with other human proteins.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant, full-length Human ERAB protein (Fusion protein).
Purification	Ammonium sulfate precipitated and dialyzed tissue culture supernatant.
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).