

Retinoic Acid Receptor gamma antibody

Cat. No. GTX102914

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse

References (1)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000

Not tested in other applications.

Calculated MW 50 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human Retinoic Acid Receptor gamma. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx102914 WB Image

Retinoic Acid Receptor gamma antibody detects RARG protein by Western blot analysis.

A. 30 µg Neuro2A whole cell lysate/extract


10 % SDS-PAGE

Retinoic Acid Receptor gamma antibody (GTx102914) dilution: 1:1000



GTx102914 WB Image

Whole cell extract (30 µg) was separated by 10% SDS-PAGE, and the membrane was blotted with Retinoic Acid Receptor gamma antibody (GTx102914) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTx213110-01) was used to detect the primary antibody.



GTx102914 WB Image

Retinoic Acid Receptor gamma antibody detects RARG protein by Western blot analysis.

A. 30 µg A549 whole cell lysate/extract

B. 30 µg HCT116 whole cell lysate/extract

10 % SDS-PAGE

Retinoic Acid Receptor gamma antibody (GTx102914) dilution: 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).