

BBOX1 antibody

Cat. No. GTX103127

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse

References (1)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 45 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 10% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.76 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human BBOX1. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx103127 WB Image

Sample (30 µg of whole cell lysate)

A: MOLT4 (GTx27912)

7.5% SDS PAGE

GTx103127 diluted at 1:500




GTx103127 IHC-P Image

BBOX1 antibody detects BBOX1 protein at cytoplasm on mouse kidney by immunohistochemical analysis.

Sample: Paraffin-embedded mouse kidney.

BBOX1 antibody (GTx103127) diluted at 1:500.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min




GTx103127 WB Image

Sample (50 µg of whole cell lysate)

A: mouse liver

10% SDS PAGE

GTx103127 diluted at 1:500



GTx103127 ICC/IF Image

Immunofluorescence analysis of methanol-fixed A431, using BBOX1 (GTx103127) antibody at 1:500 dilution.



For full product information, images and publications, please visit our [website](https://www.genetex.com).