

Apolipoprotein D antibody [8CD6]

Cat. No. GTX10336

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, ELISA
Reactivity	Human

Package
100 µl

Applications

Application Note

IHC-P: Use at a dilution of 1:200 for biotin based detection systems, for optimal staining, the sections should be pretreated with an antigen unmasking solution, and incubate for 20-60 minutes at room temperature Westernblot: Use at a dilution of 1:1,000-1:10,000 Optimal dilutions/concentrations should be determined by the end user.

Calculated MW	21 kDa. (Note)
Product Note	recognizes Apolipoprotein-D, a 24 kDa progesterone binding glycoprotein

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length protein (Human)
Purification	Tissue culture supernatant TCS
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).