

## AKAP12 antibody

Cat. No. GTX10346

Host	Sheep
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IP
Reactivity	Rat

Package  
125 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	2-5 µg/ml
IP	2-5 µg/ml

**Note : Using in the presence of 100 ng/ml of human recombinant eotaxin-2.**

Not tested in other applications.

Calculated MW 181 kDa. ([Note](#))

## Properties

Form	Liquid
Buffer	PBS
Preservative	0.08% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	rat SSeCKS.
Purification	IgG fraction
Conjugation	Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).