

MCP1 / CCL2 antibody [24822]

Cat. No. GTX10390

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, Neutralizing /Inhibition
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
Neutralizing /Inhibition	Assay dependent

Not tested in other applications.

Calculated MW 11 kDa. ([Note](#))

Product Note

The antibody will neutralize the biological activity of recombinant human MCP-1. This antibody does not exhibit cross-reactivity with recombinant human RANTES, recombinant human MIP-1 α , or recombinant mouse MIP-1 α .

Properties

Form	Liquid
Buffer	PBS carbohydrates
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	recombinant human MCP-1/CCL2 expressed in E. coli.
Purification	Purified immunoglobulin
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).