

IRAK1BP1 antibody [N1C3]

Cat. No. GTX104392

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:100-1:1000


Not tested in other applications.

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 10% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human IRAK1BP1. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX104392 IHC-P Image

IRAK1BP1 antibody [N1C3] detects IRAK1BP1 protein at cytoplasm on human placenta by immunohistochemical analysis.

Sample: Paraffin-embedded placenta.

IRAK1BP1 antibody [N1C3] (GTX104392) dilution: 1:100.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX104392 IHC-P Image

IRAK1BP1 antibody [N1C3] detects IRAK1BP1 protein at cytoplasm in mouse kidney by immunohistochemical analysis.

Sample: Paraffin-embedded mouse kidney.

IRAK1BP1 antibody [N1C3] (GTX104392) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTX104392 IHC-P Image

IRAK1BP1 antibody [N1C3] detects IRAK1BP1 protein at cytoplasm in rat intestine by immunohistochemical analysis.

Sample: Paraffin-embedded rat intestine.

IRAK1BP1 antibody [N1C3] (GTX104392) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](#).