

Bif-1 antibody [C2C3], C-term

Cat. No. GTX104636

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 41 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 10% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.57 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human Bif-1. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX104636 WB Image

Sample (30 ug of whole cell lysate)


A: NIH-3T3

B: JC

C: BCL-1

10% SDS PAGE

GTX104636 diluted at 1:1000



GTX104636 ICC/IF Image


Bif-1 antibody [C2C3], C-term detects Bif-1 protein at cytoplasm by immunofluorescent analysis.

Sample: HeLa cells were fixed in 4% paraformaldehyde at RT for 10 min.

Green: Bif-1 protein stained by Bif-1 antibody [C2C3], C-term (GTX104636) diluted at 1:200.

Blue: Hoechst 33342 staining.

Scale bar = 10 μm




GTX104636 IHC-P Image

Bif-1 antibody detects Bif-1 protein at cytosol on human ovarian carcinoma by immunohistochemical analysis.

Sample: Paraffin-embedded human ovarian carcinoma.

Bif-1 antibody (GTX104636) dilution: 1:500.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



GTX104636 WB Image

Sample(30 μg of whole cell lysate)

A:H1299

7.5% SDS PAGE

GTX104636 diluted at 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).