

RPL37 antibody

Cat. No. GTX104688

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000

Not tested in other applications.

Calculated MW 11 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Carrier-protein conjugated synthetic peptide encompassing a sequence within the center region of human RPL37. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX104688 WB Image

RPL37 antibody detects RPL37 protein by western blot analysis.


A. 30 µg Jurkat whole cell lysate/extract

B. 30 µg Raji whole cell lysate/extract

C. 30 µg NCI-H929 whole cell lysate/extract

15 % SDS-PAGE

RPL37 antibody (GTX104688) dilution: 1:1000



GTX104688 ICC/IF Image

RPL37 antibody detects RPL37 protein at cytoplasm by immunofluorescent analysis.

Sample: HeLa cells were fixed in 4% paraformaldehyde at RT for 15 min.

Green: RPL37 protein stained by RPL37 antibody (GTX104688) diluted at 1:500.

Blue: Hoechst 33342 staining.


GTX104688 WB Image

RPL37 antibody detects RPL37 protein by western blot analysis.

A. 30 µg 293T whole cell lysate/extract

B. 30 µg HeLa whole cell lysate/extract

15 % SDS-PAGE

RPL37 antibody (GTX104688) dilution: 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).