

IRF6 antibody [C3], C-term

Cat. No. GTX104862

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Rat

References (3)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:50-1:500
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 53 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 10% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human IRF6. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTXT104862 WB Image


Western Blot analysis of IRF6 expression in transfected 293T cell line by IRF6 polyclonal antibody.

A: Non-transfected lysate.

B: IRF6 transfected lysate.


12% SDS PAGE

GTXT104862 diluted at 1:500



GTXT104862 ICC/IF Image

Immunofluorescence analysis of methanol-fixed mouse ESC, using IRF6(GTXT104862) antibody at 1:100 dilution.



GTXT104862 IHC-P Image

IRF6 antibody [C3], C-term detects IRF6 protein at cytoplasm and nucleus in rat liver by immunohistochemical analysis.

Sample: Paraffin-embedded rat liver.

IRF6 antibody [C3], C-term (GTXT104862) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](https://www.genetex.com).