

CTDSP2 antibody

Cat. No. GTX104889

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse, Rat

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
Not tested in other applications.	

Calculated MW 31 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human CTDSP2. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX104889 WB Image

Sample (50 ug of whole cell lysate)

A: Mouse brain

12% SDS PAGE

GTX104889 diluted at 1:1000




GTX104889 WB Image

CTDSP2 antibody detects CTDSP2 protein by Western blot analysis.

A. 30 µg PC-12 whole cell lysate/extract

12 % SDS-PAGE

CTDSP2 antibody (GTX104889) dilution: 1:1000



GTX104889 WB Image

Sample (30 ug of whole cell lysate)

A: 293T

12% SDS PAGE

CTDSP2 antibody

GTX104889 diluted at 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).