

## LTBR antibody

Cat. No. GTX10493

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, FCM
Reactivity	Human

## Package

50 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.1-0.2 µg/ml
FCM	2.5 µg

**Note : using 1ng/lane rhTNF RII/TNFRSF1B Fc Chimera under non-reducing conditions**

Not tested in other applications.

**Calculated MW** 47 kDa. ( [Note](#) )

## Properties

Form	Liquid
Buffer	PBS carbohydrates
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	recombinant human lymphotoxin β receptor extracellular domain expressed in NS0 cells.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).