

TAF7 antibody

Cat. No. GTX105174

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IP
Reactivity	Human

★ ★ ★ ★ ★ Review (1)

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IP	1:100-1:500

Not tested in other applications.

Calculated MW 40 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.8 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human TAF7. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx105174 IP Image

TAF7 antibody immunoprecipitates TAF7 protein in IP experiments. IP Sample: 293T whole cell lysate/extract A. 40 µg 293T whole cell lysate/extract B. Control with 2 µg of preimmune rabbit IgG C. Immunoprecipitation of TAF7 protein by 2 µg of TAF7 antibody (GTx105174) 7.5% SDS-PAGE The immunoprecipitated TAF7 protein was detected by TAF7 antibody (GTx105174) diluted at 1:1000. EasyBlot anti-rabbit IgG (GTx221666-01) was used as a secondary reagent.


GTx105174 WB Image

TAF7 antibody detects TAF7 protein by western blot analysis.
A. 30 µg HeLa whole cell lysate/extract
B. 30 µg HepG2 whole cell lysate/extract
10 % SDS-PAGE
TAF7 antibody (GTx105174) dilution: 1:1000


GTx105174 WB Image

Sample (30 ug of whole cell lysate)
A: A431
B: Jurkat
10% SDS PAGE
GTx105174 diluted at 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).