

FADD antibody [FD19]

Cat. No. GTX10519

Host	Mouse
Clonality	Monoclonal
Isotype	lgG1
Applications	WB, IP, ELISA
Reactivity	Human

Package 100 μΙ

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.5-1 μg/ml
IP	Assay dependent
ELISA	Assay dependent
Not tested in other applications.	

Calculated MW 23 kDa. (Note)

Properties	
Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	~2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	purified recombinant human FADD.
Purification	Purified immunoglobulin
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our website.

Date 2025 / 11 / 07 Page 1 of 1

€ 886-3-6208988 886-3-6208989 infoasia@genetex.com