

MMP17 antibody [N2C2], Internal

Cat. No. GTX105268

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human

 Review (1)Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1500

Not tested in other applications.

Calculated MW 67 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.74 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human MMP17. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX105268 IHC-P Image

MMP17 antibody [N2C2], Internal detects MMP17 protein at cytosol and membrane on human gastric cancer by immunohistochemical analysis.

Sample: Paraffin-embedded gastric cancer.


MMP17 antibody [N2C2], Internal (GTX105268) dilution: 1:500.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



GTX105268 ICC/IF Image

Immunofluorescence analysis of methanol-fixed HeLa, using MMP17(GTX105268) antibody at 1:200 dilution.



GTX105268 IHC-P Image


MMP17 antibody [N2C2], Internal detects MMP17 protein at cytosol on human gastrica carcinoma by immunohistochemical analysis.

Sample: Paraffin-embedded human gastrica carcinoma.

MMP17 antibody [N2C2], Internal (GTX105268) dilution: 1:1500.

Scale bar = 100 μm.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



GTX105268 WB Image

Sample (30 ug of whole cell lysate)

A: Hela

7.5% SDS PAGE

GTX105268 diluted at 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).