

## RANTES antibody

**Cat. No. GTX10590**

<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG
<b>Applications</b>	WB, IHC-Fr, Neutralizing /Inhibition
<b>Reactivity</b>	Human

Package

25 µg

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.1 µg/ml
IHC-Fr	5-15 µg/ml
Neutralizing /Inhibition	0.1-0.4 µg/ml

Not tested in other applications.

**Calculated MW** 10 kDa. ([Note](#))**Product Note**

Based on ELISA, the antibody shows ~40% cross-reactivity with recombinant mouse RANTES, ~20% cross-reactivity with recombinant cotton rat RANTES, and ~1% cross-reactivity with recombinant feline RANTES.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS, 5% Trehalose
<b>Preservative</b>	No preservatives
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	Batch dependent (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	recombinant human RANTES expressed in E. coli.
<b>Purification</b>	Purified by affinity chromatography
<b>Conjugation</b>	Unconjugated



For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

**Note**

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).