

Lymphotactin antibody

Cat. No. GTX10607

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA, Neutralizing/Inhibition
Reactivity	Human

Package 25 μg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution	
WB	1-2 μg/ml	
ELISA	0.5-1.0 μg/mL	
Neutralizing/Inhibition	Assay dependent	
Not tested in other applications.		
Calculated MW	13 kDa. (<u>Note</u>)	
Product Note	The antibody will neutralize the biological activity of human lymphotactin.	
Properties		
Form	Liquid	
Buffer	PBS	
Preservative	No preservatives	
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.	
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)	
Immunogen	purified recombinant human lymphotactin expressed in Escherichia coli.	
Purification	Purified by affinity chromatography	



Conjugation

For full product information, images and publications, please visit our <u>website</u>.

Unconjugated

Date 2025 / 11 / 04 Page 1 of 2



For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 04 Page 2 of 2