

ZDHHC13 antibody [C3], C-term

Cat. No. GTX106273

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000

Not tested in other applications.

Calculated MW 71 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 10% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human ZDHHC13. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX106273 WB Image

Sample(30 µg of whole cell lysate)
A:HeLa S3 (GTX14654)
10% SDS PAGE
GTX106273 diluted at 1:500



GTX106273 ICC/IF Image

ZDHHC13 antibody [C3], C-term detects ZDHHC13 protein at cytoplasm by immunofluorescent analysis.
Sample: HeLa cells were fixed in -20°C 100% MeOH for 5 min.
Green: ZDHHC13 protein stained by ZDHHC13 antibody [C3], C-term (GTX106273) diluted at 1:500.
Blue: Hoechst 33343 staining.



For full product information, images and publications, please visit our [website](#).