

KPNA4 antibody [C3], C-term

Cat. No. GTX106325

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Rat

References (2)

★★★★★ Review (1)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 58 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.04 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human KPNA4. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx106325 WB Image


Sample (30 ug of whole cell lysate)

A: HeLa

B: Hep G2

7.5% SDS PAGE


GTx106325 diluted at 1:1000



GTx106325 IHC-P Image

Immunohistochemical analysis of paraffin-embedded H1299 xenograft, using KPNA4(GTx106325) antibody at 1:500 dilution.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



GTx106325 WB Image

Sample (30 ug of whole cell lysate)


A: 293T

B: HepG2(GTx27900)

C: Molt-4 (GTx27912)

7.5% SDS PAGE

GTx106325 diluted at 1:1500



GTx106325 ICC/IF Image

Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using QIP1 (GTx106325) antibody at 1:200 dilution.



For full product information, images and publications, please visit our [website](https://www.genetex.com).