

Coronin 1A antibody [C3], C-term

Cat. No. GTX106424

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, IHC-Fr
Reactivity	Human, Mouse, Rat

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000
IHC-Fr	1:100-1:1000

Not tested in other applications.

Calculated MW 51 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 10% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human Coronin 1A. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES

GTX106424 IHC-P Image

Coronin 1A antibody [C3], C-term detects Coronin 1A protein at cytoplasm in mouse duodenum by immunohistochemical analysis.

Sample: Paraffin-embedded mouse duodenum.

Coronin 1A antibody [C3], C-term (GTX106424) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX106424 IHC-P Image

Coronin 1A antibody [C3], C-term detects Coronin 1A protein at cytoplasm in rat brain by immunohistochemical analysis.

Sample: Paraffin-embedded rat brain.

Coronin 1A antibody [C3], C-term (GTX106424) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX106424 IHC-P Image

Coronin 1A antibody [C3], C-term detects Coronin 1A protein at cytoplasm in rat intestine by immunohistochemical analysis.

Sample: Paraffin-embedded rat intestine.

Coronin 1A antibody [C3], C-term (GTX106424) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTX106424 IHC-P Image

Coronin 1A antibody [C3], C-term detects Coronin 1A protein at cytoplasm in mouse brain by immunohistochemical analysis.

Sample: Paraffin-embedded mouse brain.

Coronin 1A antibody [C3], C-term (GTX106424) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](#).