

# Signal+ for Immunostaining - Solution B

**Cat. No. GTX10666**

**Applications**

ICC/IF, IHC

**Package**

20 ml

**PRODUCT**

**Summary**

Signal+ for Immunostaining - Solution B

**Applications**

**Application Note**

Signal+ for Immunostaining Solutions A and B exhibit different enhancing effects, depending on antigens and antibodies. These solutions can be used independently. However, both solutions should be explored in order to select the more suitable one for your particular assay system. A Sample Kit composed of both Solution A and B is available as GTX49998.

**Properties**

**Form**

Liquid

**Storage**

Store at 4°C.

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).