

GM-CSF antibody

Cat. No. GTX10669

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA
Reactivity	Rat

Package

25 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.1-0.2 µg/ml
ELISA	0.5-1.0 µg/ml

Note : At 0.5-2 µg/mL the antibody will block 50% of the binding of 5 ng/mL recombinant human HGF to 0.1 µg immobilized recombinant human HGF R/c-MET Fc chimera. >90% of the binding will be blocked by 10 µg/mL of the antibody.

Not tested in other applications.

Calculated MW	17 kDa. (Note)
Product Note	By immunoblotting and ELISA, the antibody shows 30% cross-reactivity with recombinant mouse GM-CSF and 10% cross-reactivity with recombinant human GM-CSF.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Purified recombinant rat GM-CSF expressed in E. coli.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 06 Page 1 of 2

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).