

TATDN1 antibody [N1C2]

Cat. No. GTX106818

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 34 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 10% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human TATDN1. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 07 Page 1 of 2

DATA IMAGES

GTX106818 IHC-P Image

TATDN1 antibody [N1C2] detects TATDN1 protein at cytosol on human ovarian carcinoma by immunohistochemical analysis.

Sample: TATDN1 antibody [N1C2] (GTX106818) dilution: 1:250.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min


GTX106818 WB Image

Sample(30 ug whole cell lysate)


A:293T

B:A431(GTX27909)

C:MOLT4 (GTX27912)

12% SDS PAGE

GTX106818 diluted at 1:1000


GTX106818 ICC/IF Image

Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using TATDN1(GTX106818) antibody at 1:200 dilution.



For full product information, images and publications, please visit our [website](https://www.genetex.com).