

U2AF1 antibody

Cat. No. GTX106854

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, IP
Reactivity	Human, Mouse

References (1)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000
IP	1:100-1:500

Not tested in other applications.

Calculated MW 28 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	0.025% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.25 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human U2AF35. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES

GTX106854 IHC-P Image

U2AF35 antibody detects U2AF35 protein at nucleus on mouse intestine by immunohistochemical analysis.

Sample: Paraffin-embedded mouse intestine.

U2AF35 antibody (GTX106854) diluted at 1:1000.


Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min


GTX106854 IP Image

Immunoprecipitation of U2AF35 protein from 293T whole cell extracts using 5 µg of U2AF35 antibody (GTX106854).

Western blot analysis was performed using U2AF35 antibody (GTX106854).

EasyBlot anti-Rabbit IgG (GTX221666-01) was used as a secondary reagent.


GTX106854 ICC/IF Image

Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using U2AF35(GTX106854) antibody at 1:200 dilution.



For full product information, images and publications, please visit our [website](https://www.genetex.com).