

CCDC68 antibody [N2C3]

Cat. No. GTX106883

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IP
Reactivity	Human

References (2)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IP	Assay dependent

Not tested in other applications.

Calculated MW 39 kDa. ([Note](#))

Product Note Knockdown/Knockout validation was supported by references (PMID:28422092).
KO/KD validation is based on published data (PMID: 28422092).

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 10% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human CCDC68. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 06 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTx106883 ICC/IF Image

Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using CCDC68(GTx106883) antibody at 1:200 dilution.



GTx106883 WB Image

Sample(30 ug whole cell lysate)
A:Raji (GTx27908)
10% SDS PAGE
GTx106883 diluted at 1:1000



GTx106883 ICC/IF Image

The data was published in the journal Nat Commun in 2017.[PMID: 28422092](https://pubmed.ncbi.nlm.nih.gov/28422092/)



For full product information, images and publications, please visit our [website](https://www.genetex.com).