

GABA A Receptor alpha 1 antibody [S95-35]

Cat. No. GTX10698

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Application	WB, ICC/IF
Reactivity	Human, Mouse, Rat

Package
100 µg

APPLICATION

Application Note

For WB: 1-10ug/ml For IF: 1-10ug/ml

Calculated MW 55 kDa. ([Note](#))

Product Note No cross-reactivity against GABA-A-R-Alpha 2 or -Alpha3.

PROPERTIES

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Fusion protein amino acids 355-394 of mouse GABA-A-R-Alpha1, accession number P62812
Purification	Protein G purified
Conjugation	Unconjugated
Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>

For full product information, images and publications, please visit our [website](#).