

GABA A Receptor alpha 1 antibody [S95-35]

Cat. No. GTX10698

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ICC/IF
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

For WB: 1-10ug/ml For IF: 1-10ug/ml

Calculated MW 55 kDa. ([Note](#))**Product Note** No cross-reactivity against GABA-A-R-Alpha 2 or -Alpha3.

Properties

Form Liquid**Buffer** PBS, 50% Glycerol**Preservative** 0.09% Sodium azide**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)**Immunogen** Fusion protein amino acids 355-394 of mouse GABA-A-R-Alpha1, accession number P62812**Purification** Protein G purified**Conjugation** Unconjugated**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).