

MPP3 antibody [N1C1]

Cat. No. GTX107608

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 66 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 10% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.39 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human MPP3. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx107608 WB Image

Sample (30 ug of whole cell lysate)
A: H1299
7.5% SDS PAGE
GTx107608 diluted at 1:500



GTx107608 IHC-P Image


Immunohistochemical analysis of paraffin-embedded SW480 xenograft, using MPP3(GTx107608) antibody at 1:500 dilution.
Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTx107608 IHC-P Image

MPP3 antibody [N1C1] detects MPP3 protein at cytoplasm and nucleus in mouse lung by immunohistochemical analysis.
Sample: Paraffin-embedded mouse lung.
MPP3 antibody [N1C1] (GTx107608) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTx107608 IHC-P Image

MPP3 antibody [N1C1] detects MPP3 protein at cytoplasm and nucleus in mouse brain by immunohistochemical analysis.
Sample: Paraffin-embedded mouse brain.
MPP3 antibody [N1C1] (GTx107608) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](https://www.genetex.com).