

## RIPK5 antibody [C1C3-2]

Cat. No. GTX107925

|              |            |
|--------------|------------|
| Host         | Rabbit     |
| Clonality    | Polyclonal |
| Isotype      | IgG        |
| Applications | WB, IHC-P  |
| Reactivity   | Human      |

Package  
100 µl, 25 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:500-1:3000         |
| IHC-P              | 1:100-1:1000         |

Not tested in other applications.

Calculated MW 105 kDa. ([Note](#))

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | 0.1M Tris, 0.1M Glycine, 10% Glycerol                                                                                                                                                                                      |
| Preservative  | 0.01% Thimerosal                                                                                                                                                                                                           |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                   |
| Immunogen     | Recombinant protein encompassing a sequence within the center region of human RIPK5. The exact sequence is proprietary.                                                                                                    |
| Purification  | Purified by antigen-affinity chromatography.                                                                                                                                                                               |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |


## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTx107925 IHC-P Image**

Immunohistochemical analysis of paraffin-embedded H1299 xenograft , using RIPK5(GTx107925) antibody at 1:100 dilution.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min


**GTx107925 WB Image**

Sample (30µg whole cell lysate)

A:Hep G2 (GTx27900)

B:MOLT4 (GTx27912)

7.5% SDS PAGE

GTx107925 diluted at 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).