

COASY antibody [N3C3-3]

Cat. No. GTX107934

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 62 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 10% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human COASY. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx107934 ICC/IF Image

Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using CoA synthase(GTx107934) antibody at 1:200 dilution.



GTx107934 IHC-P Image

CoA synthase antibody [N3C3-3] detects CoA synthase protein at cytoplasm in mouse kidney by immunohistochemical analysis.

Sample: Paraffin-embedded mouse kidney.

CoA synthase antibody [N3C3-3] (GTx107934) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTx107934 IHC-P Image

CoA synthase antibody [N3C3-3] detects CoA synthase protein at cytoplasm in mouse prostate by immunohistochemical analysis.

Sample: Paraffin-embedded mouse prostate.

CoA synthase antibody [N3C3-3] (GTx107934) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTx107934 WB Image

Sample (30µg whole cell lysate)

A:HeLa S3(GTx14654)

7.5% SDS PAGE

GTx107934 diluted at 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).