

VAM1 antibody [N1], N-term

Cat. No. GTX108010

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

References (3)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 61 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 10% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Carrier-protein conjugated synthetic peptide encompassing a sequence within the N-terminus region of human VAM1. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx108010 WB Image

Sample (50 ug of whole cell lysate)

A: Rat brain

7.5% SDS PAGE

GTx108010 diluted at 1:1000




GTx108010 IHC-P Image

VAM1 antibody [N1], N-term detects VAM1 protein at on mouse intestine by immunohistochemical analysis.

Sample: Paraffin-embedded mouse intestine.

VAM1 antibody [N1], N-term (GTx108010) dilution: 1:500.


Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



GTx108010 IHC-P Image

Immunohistochemical analysis of paraffin-embedded human colon, using VAM1(GTx108010) antibody at 1:100 dilution.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



GTx108010 WB Image

Sample (50 ug of whole cell lysate)

A: mouse brain

7.5% SDS PAGE

GTx108010 diluted at 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).