

## Ephrin B3 antibody [C2C3], C-term

Cat. No. GTX108485

|              |                   |
|--------------|-------------------|
| Host         | Rabbit            |
| Clonality    | Polyclonal        |
| Isotype      | IgG               |
| Applications | WB                |
| Reactivity   | Human, Mouse, Rat |

Package  
100 µl, 25 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:500-1:3000         |

Not tested in other applications.

Calculated MW 36 kDa. ( [Note](#) )

Product Note This antibody was raised against human Ephrin B3 Intracellular domain.

## Properties


|               |                                                                                                                                                                                                                                                                                                                                             |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                                                                                                                                      |
| Buffer        | 0.1M Tris, 0.1M Glycine, 10% Glycerol                                                                                                                                                                                                                                                                                                       |
| Preservative  | 0.01% Thimerosal                                                                                                                                                                                                                                                                                                                            |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.                                                                                                                  |
| Concentration | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                                                                                                                                    |
| Immunogen     | Synthetic peptide encompassing a sequence within the Intracellular domain of human Ephrin B3. The exact sequence is proprietary.                                                                                                                                                                                                            |
| Purification  | Purified by antigen-affinity chromatography.                                                                                                                                                                                                                                                                                                |
| Conjugation   | Unconjugated                                                                                                                                                                                                                                                                                                                                |
| Note          | <p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p> |



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 07 Page 1 of 2

## DATA IMAGES




### GTX108485 WB Image

Rat tissue extract (50  $\mu$ g) was separated by 12% SDS-PAGE, and the membrane was blotted with Ephrin B3 antibody [C2C3], C-term (GTX108485) diluted at 1:500. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.



### GTX108485 WB Image

Various whole cell extracts (30  $\mu$ g) were separated by 12% SDS-PAGE, and the membrane was blotted with Ephrin B3 antibody [C2C3], C-term (GTX108485) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primar



### GTX108485 WB Image

Mouse tissue extract (50  $\mu$ g) was separated by 12% SDS-PAGE, and the membrane was blotted with Ephrin B3 antibody (GTX108485) diluted at 1:500. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).