

TRAM1 antibody [C2C3], C-term

Cat. No. GTX108535

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 43 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 10% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human TRAM1. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTx108535 WB Image

Sample (30 ug of whole cell lysate)

A: H1299

10% SDS PAGE

GTx108535 diluted at 1:1000



GTx108535 IHC-P Image

TRAM1 antibody [C2C3], C-term detects TRAM1 protein at cytoplasm on AGS xenograft by immunohistochemical analysis.

Sample: Paraffin-embedded AGS xenograft.

TRAM1 antibody [C2C3], C-term (GTx108535) dilution: 1:500.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



GTx108535 IHC-P Image

TRAM1 antibody [C2C3], C-term detects TRAM1 protein at cytoplasm on mouse intestine by immunohistochemical analysis.

Sample: Paraffin-embedded mouse intestine.

TRAM1 antibody [C2C3], C-term (GTx108535) diluted at 1:500.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min


GTx108535 ICC/IF Image

Immunofluorescence analysis of methanol-fixed HeLa, using TRAM1(GTx108535) antibody at 1:500 dilution.



For full product information, images and publications, please visit our [website](https://www.genetex.com).