

Human Rheb protein, His tag

Cat. No. GTX108597-pro

Species Human

Package

100 µg, 25 µg

PROPERTIES

Form Liquid

Buffer PBS, 0.1% SDS

Preservative No Preservative

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration Batch dependent (Please refer to the vial label for the specific concentration.)

Region/Sequence N-terminal His tagged partial fragment of NP_005605


Expression System E. coli

Purity >95%

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

DATA IMAGES



GTX108597-pro Image

10µg of GTX108597-pro RHEB recombinant protein analyzed using SDS-PAGE and stained with coomassie blue and captured by black and white camera.



For full product information, images and publications, please visit our [website](https://www.genetex.com).