

PPP6C antibody

Cat. No. GTX108968

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, IHC-Wm
Reactivity	Human, Zebrafish

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	1:100-1:1000
IHC-Wm	Assay dependent

Not tested in other applications.

Calculated MW 35 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 10% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.51 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length human PPP6C Recombinant protein.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX108968 WB Image

Sample (30 µg of whole cell lysate)
A: Zebrafish brain
10% SDS PAGE
GTX108968 diluted at 1:1000




GTX108968 IHC-Wm Image

Immunohistochemical analysis (whole mount) of zebrafish embryo, using PPP6C antibody (GTX108968) at 1:200 dilution.



GTX108968 WB Image

Sample (30 ug of whole cell lysate)
A: HeLa
10% SDS PAGE
GTX108968 diluted at 1:1000



GTX108968 IHC-P Image

Immunohistochemical analysis of paraffin-embedded PC14 xenograft, using PPP6C(GTX108968) antibody at 1:500 dilution.
Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](https://www.genetex.com).