

Human L-FABP protein, His tag

Cat. No. GTX109142-pro

Species	Human
---------	-------


Package

100 µg, 25 µg

Properties

Form	Liquid
Buffer	PBS, 0.1% SDS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Region/Sequence	N-terminal His tagged full length protein of NP_001434
Expression System	E. coli
Purity	>95%
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

DATA IMAGES



GTX109142-pro Image

10µg of GTX109142-pro FABP1 (liver) recombinant protein analyzed using SDS-PAGE and stained with coomassie blue and captured by black and white camera.



For full product information, images and publications, please visit our [website](#).