

Human C/EBP epsilon protein, His tag

Cat. No. GTX109155-pro

Species Human

Package

100 µg, 25 µg

Properties

Form Liquid

Buffer PBS, 0.1% SDS

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration Batch dependent (Please refer to the vial label for the specific concentration.)

Region/Sequence N-terminal His tagged full length protein of NP_001796

Expression System E. coli

Purity >90%

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 08 Page 1 of 1