

LHX8 antibody

Cat. No. GTX109165

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IP
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000-1:10000
IP	1:100-1:500

Not tested in other applications.

Calculated MW 39 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 10% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length human LHX8 Recombinant protein.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX109165 IP Image

LIM homeobox protein 8 antibody immunoprecipitates LIM homeobox protein 8 protein in IP experiments.
 IP Sample: 293T whole cell lysate/extract A. 40 µg 293T whole cell lysate/extract B. Control with 2 µg of preimmune rabbit IgG C. Immunoprecipitation of LIM homeobox protein 8 protein by 2 µg of LIM homeobox protein 8 antibody (GTX109165) 10% SDS-PAGE The immunoprecipitated LIM homeobox protein 8 protein was detected by LIM homeobox protein 8 antibody (GTX109165) diluted at 1:1000. EasyBlot anti-rabbit IgG (GTX221666-01) was used as a secondary reagent.



GTX109165 WB Image

Sample (30 ug of whole cell lysate)
 A: Molt-4 (GTX27912)
 10% SDS PAGE
 GTX109165 diluted at 1:5000



For full product information, images and publications, please visit our [website](https://www.genetex.com).