

Telethonin antibody [N1C3]

Cat. No. GTX109265

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 19 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 10% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length human Telethonin Recombinant protein.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX109265 IHC-P Image

Telethonin antibody [N1C3] detects Telethonin protein at cytoplasm in Huh7 xenograft by immunohistochemical analysis.

Sample: Paraffin-embedded Huh7 xenograft.

Telethonin antibody [N1C3] (GTX109265) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX109265 WB Image

Telethonin antibody detects Telethonin protein by western blot analysis. Mouse tissue extracts (50 µg) was separated by 15% SDS-PAGE, and the membrane was blotted with Telethonin antibody (GTX109265) diluted at 1:500.



For full product information, images and publications, please visit our [website](https://www.genetex.com).