

DUSP26 antibody [N1C1]

Cat. No. GTX109283

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IP
Reactivity	Human, Mouse, Rat

References (6)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IP	Assay dependent

Not tested in other applications.

Calculated MW	24 kDa. (Note)
Product Note	KO/KD validation is based on published data (PMID: 24548998).

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	0.025% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.52 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Carrier-protein conjugated synthetic peptide encompassing a sequence within the center region of human DUSP26. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 06 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTx109283 ICC/IF Image

Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using DUSP26(GTx109283) antibody at 1:200 dilution.



GTx109283 WB Image

Sample (30 ug of whole cell lysate)
A: IMR32
12% SDS PAGE
GTx109283 diluted at 1:500



GTx109283 WB Image

Sample (50 ug of whole cell lysate)
A: mouse brain
12% SDS PAGE
GTx109283 diluted at 1:500



For full product information, images and publications, please visit our [website](https://www.genetex.com).