

SmarcAL1 antibody

Cat. No. GTX109468

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IP
Reactivity	Human, Mouse, Rat

References (2)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IP	1:100-1:500

Not tested in other applications.

Calculated MW 106 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	0.025% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.07 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human SmarcAL1. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx109468 ICC/IF Image

Immunofluorescence analysis of paraformaldehyde-fixed A549, using SmarAL1(GTx109468) antibody at 1:500 dilution.


GTx109468 IP Image

SmarAL1 antibody immunoprecipitates SMARAL1 protein in IP experiments.

IP samples: A549 whole cell extract


A. Control with 4 µg of preimmune Rabbit IgG

B. Immunoprecipitation of SMARAL1 protein by 4 µg SmarAL1 antibody (GTx109468)


5 % SDS-PAGE

The immunoprecipitated SMARAL1 protein was detected by SmarAL1 antibody (GTx109468) diluted at 1:1000.

[EasyBlot anti-rabbit IgG (GTx221666-01) was used as a secondary reagent]


GTx109468 WB Image

Various whole cell extracts (30 µg) were separated by 5% SDS-PAGE, and the membrane was blotted with SmarAL1 antibody (GTx109468) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTx213110-01) was used to detect the primary antibody.


GTx109468 WB Image

Mouse tissue extract (50 µg) was separated by 7.5% SDS-PAGE, and the membrane was blotted with SmarAL1 antibody (GTx109468) diluted at 1:500. The HRP-conjugated anti-rabbit IgG antibody (GTx213110-01) was used to detect the primary antibody.



For full product information, images and publications, please visit our [website](http://www.genetex.com).