

DLL4 antibody [N1N2], N-term

Cat. No. GTX109649

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human

References (2)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000

Not tested in other applications.

Calculated MW 75 kDa. ([Note](#))


Product Note This antibody was raised against human DLL4 Extracellular domain.

Properties


Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 10% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide encompassing a sequence within the Extracellular domain of human DLL4. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX109649 WB Image


Whole cell extract (30 µg) was separated by 7.5% SDS-PAGE, and the membrane was blotted with DLL4 antibody [N1N2], N-term (GTX109649) diluted at 1:2000.


GTX109649 WB Image

Sample (30 ug of whole cell lysate)
A: Raji
7.5% SDS PAGE
DLL4 antibody
GTX109649 diluted at 1:2000


GTX109649 WB Image

The data was published in the Biol Res in 2020. [PMID: 32046779](https://pubmed.ncbi.nlm.nih.gov/32046779/)


GTX109649 WB Image

The data was published in the Biol Res in 2020. [PMID: 32046779](https://pubmed.ncbi.nlm.nih.gov/32046779/)



For full product information, images and publications, please visit our [website](https://www.genetex.com).