

PLGF2 antibody [62526]

Cat. No. GTX10967

Host	Rat
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ELISA, Neutralizing /Inhibition
Reactivity	Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1-2 µg/ml
ELISA	0.5-1.0 µg/ml
Neutralizing /Inhibition	Assay dependent

Note : using peptide coating at 0.5-1 µg/well in 0.05 M carbonate/bicarbonate buffer, pH 9.6, and pNPP as substrate

Not tested in other applications.

Calculated MW 18 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 5% Trehalose
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	purified recombinant mouse placenta growth factor 2 expressed in insect Sf 21 cells.
Purification	Purified immunoglobulin
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).



For full product information, images and publications, please visit our [website](#).