

# Corticosterone antibody

**Cat. No. GTX10996**

<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG
<b>Applications</b>	RIA
<b>Reactivity</b>	Species independent

**Package**  
100 test

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
RIA	Assay dependent

Not tested in other applications.

## Properties

<b>Form</b>	Lyophilized powder
<b>Buffer</b>	Reconstitute with 5ml buffer. Buffer contents: 0.05M Tris-HCl pH8.0, 0.1M Sodium Chloride, 0.1% BSA
<b>Preservative</b>	0.1% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
<b>Immunogen</b>	Corticosterone-21-Thyroglobulin
<b>Purification</b>	Unpurified
<b>Conjugation</b>	Unconjugated

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).