

Corticosterone antibody

Cat. No. GTX10996

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	RIA
Reactivity	Species independent

Package
100 test

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
RIA	Assay dependent

Not tested in other applications.

Properties

Form	Lyophilized powder
Buffer	Reconstitute with 5ml buffer. Buffer contents: 0.05M Tris-HCl pH8.0, 0.1M Sodium Chloride, 0.1% BSA
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Immunogen	Corticosterone-21-Thyroglobulin
Purification	Unpurified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).