

Hepatic lipase antibody [N1C1]

Cat. No. GTX110336

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 56 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 10% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human Hepatic lipase. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx110336 IHC-P Image

Hepatic lipase antibody [N1C1] detects Hepatic lipase protein at cytoplasm in Hep G2 xenograft by immunohistochemical analysis.

Sample: Paraffin-embedded Hep G2 xenograft.

Hepatic lipase antibody [N1C1] (GTx110336) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTx110336 WB Image

Hepatic lipase antibody [N1C1] detects Hepatic lipase [N1C1] protein by western blot analysis. Whole cell extracts (30 µg) was separated by 10% SDS-PAGE, and the membrane was blotted with Hepatic lipase antibody [N1C1] (GTx110336) diluted at 1:500.



For full product information, images and publications, please visit our [website](https://www.genetex.com).